

**The PRECISION you
require with the
VERSATILITY you need
to be successful**



Radiofrequency Energy Sources & Accessory Products

Dermatology

Ear, Nose, and Throat (ENT)

Facial Plastic

General Surgery

Gynecology

Ophthalmology / Oculoplastics

Oral / Maxillofacial Surgery

Plastic Reconstructive

Primary Care

Podiatry

Proctology




Radiofrequency (CRF) Energy Sources

Five Distinct Waveforms	4
Pellevé®S5	5
Surgitron® Dual RF™ S5	6
Surgitron®Dual RF™ and Surgitron® Dual/EMC	7
AcuSect™	8
Surgitron®FFPF EMC	8




RF Accessories

Neutral Plates	9
Footswitches	10
Cables	11
Surg-e®-Doc™	12
EMC Device Accessories	13
Handpieces	15
Bipolar Forceps	17
Bayonet Bipolar Forceps	19
ACe-Tip™ Non-Stick Forceps	20



RF Electrodes

Sterile/Single-Use Electrodes	23
Empire® microIncision™ Needles	24
Reusable Electrodes	26
ACe-Tip™ Electrodes	31
Electrode Sets	32
Electrode Accessories	34



Educational Materials

Patient Education Brochures	35
Pellevé® Marketing Kits	35
Laser Office Pull Up Banners	35
Pellevé® Posters	35
Pellevé® Before & After Book	35



Smoke Evacuation System

Accessories	37
-------------	----

Multi Specialty Accessories

Otorhinolaryngology (ENT)	39
General Surgery/Proctology	46
Gynecology	47
Plastic Surgery	49
Dermatology	51
Pellevé® Wrinkle Reduction System	52
Ophthalmology/Oculoplastics	53
Family Practice/Podiatry	55

Index

4

4
5
6
7
8
8

9

9
10
11
12
13
15
17
19
20

22

23
24
26
31
32
34

35

35
35
35
35
35

36

37

38

39
46
47
49
51
52
53
55

56

For orders and/or pricing information, please call ellman® at: (800) 835-5355
(516) 594-3333 • fax (516) 279-4483 • www.ellman.com

RADIOFREQUENCY ENERGY SOURCES

5 Distinct Waveforms

1

Cut (Fully Filtered)




- Micro-smooth cutting
- Negligible lateral heat
- Minimal cellular destruction
- Best cosmetic results. Fastest healing^{1,2}.
- Ideal for skin incision and biopsy

ellman's radiofrequency technology has received worldwide acclaim and recognition from leading surgeons, clinics, universities and authors.

2

Blend (Cut/Coag) (Fully Rectified)



- Cutting with hemostasis
- Ideal for sub-cutaneous tissue dissection and planing. Especially useful in vascular areas while producing minimal amounts of lateral heat and tissue damage

3

Coag (Hemo) (Partially Rectified)



- Coagulation / Shrinkage
- Hemostasis with controlled penetration
- Ideal for cutting with hemostatic control

4

Fulguration




- Maximum penetration and hemostasis
- Ideal for intentional tissue destruction

5

Bipolar


- Pinpoint, micro-coagulation
- Minimal tissue adherence to forceps
- Minimal charring or tissue necrosis
- Ideal for coagulation in and around critical anatomy



Clinical Citations

- 1.Bridenstine, J.B., Derm Surgery (1998); vol 24, p397-400
- 2.Aferzon, M. Derm Surgery (2002); vol 28, p735-738

TRANSMISSION ELECTRON MICROSCOPY SHOWING LATERAL THERMAL DAMAGE TO HUMAN OVIDUCTS BY DIFFERENT INSTRUMENTS



Reference - Olivar, AC, et al, Ann Clin Lab Sci. 1999 Oct-Dec; 29(4): p281-5.

The ellman® radiofrequency RF source represents advanced radiowave technology that provides outstanding surgical control, precision and versatility. Clinical benefits include reduced post-operative discomfort, minimal scar tissue formation, enhanced healing and excellent cosmetic results.

For orders and/or pricing information, please call ellman® at: (800) 835-5355
(516) 594-3333 • fax (516) 279-4483 • www.ellman.com

RADIOFREQUENCY ENERGY SOURCES

Pellevé® S5



RF Energy Source

Versatility with the proven Surgitron® platform

The Pellevé® S5 System is a newly enhanced Pellevé-branded version of our proven Surgitron® Dual RF™ S5, providing maximum control in precision cutting and energy delivery, plus more versatility than other energy-based technologies, including:

- GlideSafe® Wrinkle Treatment Handpieces for skin tightening (see page 52)
- Hundreds of surgical electrodes and accessories for most specialties
- Ideal for office and clinic environments

Powered by advanced radiofrequency technology, the patented Pellevé® S5 energy source provides unparalleled surgical PRECISION in soft tissue cutting and coagulation.

The high frequency of 4.0 MHz minimizes thermal tissue damage along the incision path.

The S5 offers complete VERSATILITY with more options for use in surgical and aesthetic procedures than other energy-based platforms and includes the Pellevé® Wrinkle Reduction System, a revolutionary tool for the non-ablative treatment of mild to moderate facial wrinkles.

The patented Pellevé® S5 redefines your surgical results with two distinct frequencies – 4.0 MHz (Monopolar) and 1.7 MHz (Bipolar) – for outstanding precision and control. The Pellevé® S5 and Surgition® Dual RF™ S5 are the only Ellman energy sources compatible with GlideSafe® Wrinkle Treatment Handpieces.

See page 36 for Cart and Accessories
See page 52 for GlideSafe® Handpieces

For orders and/or pricing information, please call ellman® at: (800) 835-5355
(516) 594-3333 • fax (516) 279-4483 • www.ellman.com

Pellevé® S5 – IEC

*Item # IEC5P-ST

Shipped complete with:

*IEC-PC110 (for 110v) power cord

*IEC-PC220 (for 220v) power cord

DF-FSC# Dual Footswitch & Cable

IEC-HP1 IEC Foot-Controlled Handpiece

IEC-3FHPB IEC 3-Button Fingerswitch Handpiece

IEC-3FHPB/D1 Sterile/Single-Use 3-Button Fingerswitch Handpiece

IEC-BC/I 3-Bipolar Sterile/Single-Use Cables

IEC-NPD IEC Neutral Plate Sterile/Single-Use

DSEP40 40 Medical Electrodes Assortment Package Sterile/Single-Use

H79 Handpiece Clip

"J" Series Bipolar Forceps (Item number determined by selection) see page 17 & 18

DIM-92-55 Instruction Manual (USA)

DIM-92-56 Instruction Manual (International)

NOTE: Included accessories may vary by country and/or by special packages.



RADIOFREQUENCY ENERGY SOURCES

Surgitron® Dual RF™ S5



RF Energy Source

Powered by advanced radiowave technology, the patented Surgitron® Dual RF™ S5 energy source provides unparalleled surgical PRECISION in soft tissue cutting and coagulation. The high frequency of 4.0 MHz minimizes thermal tissue damage along the incision path¹.

The S5 offers complete VERSATILITY, with more options for use in surgical and aesthetic procedures than other energy-based platforms,² and including Pellevé™ accessories, revolutionary tools for the non-ablative treatment of mild to moderate facial wrinkles.

The patented Surgitron® Dual RF™ S5 redefines your surgical results with two distinct frequencies – 4.0 MHz (Monopolar) and 1.7 MHz (Bipolar) – for outstanding precision and control.

The Surgitron® Dual RF™ S5 and the Pelleve® S5 RF energy sources are the only Ellman energy sources compatible with GlideSafe® Wrinkle Treatment Handpieces.

Surgitron® Dual RF™ S5 – IEC

*Item # IEC5-S30

Shipped complete with:

*IEC-PC110 (for 110v) power cord

*IEC-PC220 (for 220v) power cord

DF-FSC#	Dual Footswitch & Cable
IEC-HP1	IEC Foot-Controlled Handpiece
IEC-3FHPB	IEC 3-Button Fingerswitch Handpiece
IEC-3FHPB/D1	Sterile/Single-Use 3-Button Fingerswitch Handpiece
IEC-BC/I	3-Bipolar Sterile/Single-Use Cables
IEC-NPD	IEC Neutral Plate Sterile/Single-Use
DSEP40	40 Medical Electrodes Assortment Package Sterile/Single-Use
H79	Handpiece Clip
"J" Series	Bipolar Forceps (Item number determined by selection) see page 17 & 18
DIM-50-48	Instruction Manual (USA)
DIM-92-52	Instruction Manual (International)

NOTE: Included accessories may vary by country and/or by special packages.

See page 36 for Cart and Accessories
See page 52 for GlideSafe® Handpieces

For orders and/or pricing information, please call ellman® at: (800) 835-5355
(516) 594-3333 • fax (516) 279-4483 • www.ellman.com

RF Energy Sources
GlideSafe® Handpieces
Pelleve® S5 RF Energy Sources

RADIOFREQUENCY ENERGY SOURCES

Surgitron® Dual RF™ 120



□□□□□□□□ ! "#

RF Energy Source 120

The patented Surgitron® Dual RF™ generator for those doctors wanting to perform advanced surgical procedures, but are not interested in wrinkle treatment. The high frequency of 4.0 MHz minimizes heat dissipation and thus cellular alteration. The system generates 120 watts of power and operates at 4.0 MHz in monopolar mode and 1.7 MHz in bipolar mode.

Surgitron® Dual/EMC 90



□□□□□□□□ ! "#

RF Energy Source 90

The system generates 90 watts and operates at 4.0 MHz in monopolar mode and 1.7 MHz in bipolar mode.

Surgitron® Dual RF™/120 IEC

*Item # IEC3A-S30

Shipped complete with:

*IEC-PC110 (for 110v) power cord

*IEC-PC220 (for 220v) power cord

DF-FSC# Dual Footswitch & Cable

IEC-HP1 IEC Foot-Controlled Handpiece

IEC-3FHPB IEC 3-Button Fingerswitch Handpiece

IEC-3FHPBD1 Sterile/Single-Use 3-Button Fingerswitch Handpiece

IEC-BC/I 3-Bipolar Sterile/Single-Use

IEC-NPD IEC Neutral Plate, Sterile/Single-Use Cables

IEC-NPC IEC Neutral Plate

DSEP40 40 Medical Electrodes Assortment Package Sterile/Single-Use

H79 Handpiece Clip

"J" Series Bipolar Forceps (Item number determined by selection) see page 17 & 18

DIM-92-51 Instruction Manual (USA)

DIM-92-52 Instruction Manual (International)

NOTE: Included accessories may vary by country and/or by special packages.

Surgitron® Dual/EMC/90

*Item # IEC3a-S30-90

Shipped complete with:
with same accessories as above.

□□□□□□□□□□□□□□□□□□□□□□

□□□□□□□□□□□□

□□□□□□□□□□□□□□□□ ! "#

For orders and/or pricing information, please call ellman® at: (800) 835-5355
(516) 594-3333 • fax (516) 279-4483 • www.ellman.com

RADIOFREQUENCY ENERGY SOURCES

AcuSect®



RF Energy Source

The AcuSect® radiosurgical energy source delivers advanced radiowave technology providing outstanding surgical control, precision and versatility. The high frequency of 4.0 MHz minimizes heat dissipation, and thus cellular alteration, while cutting and coagulating soft tissues. AcuSect® is the ideal solution for your office-based minor surgical procedures. The system generates 50 watts of power and operates at 4.0 MHz in 3 monopolar modes (Cut/Cut-Coag/Hemo) and 1 bipolar mode (Hemo).

Surgitron® FFPF EMC



RF Energy Source

The ellman® Surgitron® FFPF EMC cuts and coagulates soft tissue using monopolar energy. This system generates 90 watts of power and operates at 3.8 MHz.

AcuSect®

*Item # ACS-S50/110

Shipped complete with:

FTSW-2PC	Footswitch 2 Pedal Compact
IEC-1FHPA	One-Button Fingerswitch HP
IEC-HP1	Foot-Controlled Handpiece
IEC-NPC	Neutral Plate
IEC-BC/1	Minor Surgery Bipolar Cable (3pcs)
"J" Series	Bipolar Forceps (Item number determined by selection) see page 17 & 18
DSEP40	Medical Electrode Assortment

NOTE: Included accessories may vary by country and/or by special packages.

Surgitron® FFPF EMC S10B

*Item # EMC 110A (for 110v) Power Cord

*Item # EMC 220A (for 220v) Power Cord

Shipped complete with:

HP1	Foot-Controlled Handpiece with Locking Chuck
H4	EMC Neutral Plate
*PC	EMC Power Cord
	Standard electrode set, reusable (7pcs)
EMCL1	Operating Manual
H79	Handpiece Clip

NOTE: Included accessories may vary by country and/or by special packages.

For orders and/or pricing information, please call ellman® at: (800) 835-5355
(516) 594-3333 • fax (516) 279-4483 • www.ellman.com

RF ACCESSORIES

Neutral Plates

IEC Neutral Plate

IEC-NPC

Reusable, coated

For Dual RF™, AcuSect™ and RadioLase® Devices
(1 pc)



With the **ellman®**

Reusable Neutral Plate,
no skin contact required.

IEC Disposable Neutral Plate-Standard Size

IEC-NPD

Single-Use

For Dual RF™, S5, AcuSect™ and RadioLase® Devices
(25 pcs)



IEC Disposable Neutral Plate-Small Size

IEC-NPD-S

Single-Use

For Dual RF™, S5, AcuSect™ and RadioLase® Devices
(25 pcs)




The **ellman®** Single-Use Neutral Plate is a sensing, pre-gelled, self-adhering, dispersive pad. The effective electrodes are two aluminum oil conductors, covered by a hydrogel conductive adhesive. The conductive hydrogel provides a low-current density return path from the patient to the RF Device while the two foils provide a means of sensing contact quality.

EMC Neutral Plate

H4

Reusable, coated

For the EMC Device
(1 pc)



For orders and/or pricing information, please call **ellman®** at: (800) 835-5355
(516) 594-3333 • fax (516) 279-4483 • www.ellman.com

RF ACCESSORIES

Footswitches

Dual Frequency Footswitch and Cable

DF-FSC#

For Surgitron® Dual RF™ and S5 Devices
(1 pc)



Single-Action Bipolar Footswitch and Cable

DF-FSC/3

For Surgitron® Dual RF™ and S5 Devices
(1 pc)



Footswitch, 2 Pedal, Compact

FTSW-2PC

For AcuSect™ and RadioLase® Devices
(1 pc)



Pneumatic Footswitch

SVFSC

For use with Surg-e®-Vac™
Smoke Evacuation System
(1 pc)



For orders and/or pricing information, please call ellman® at: (800) 835-5355
(516) 594-3333 • fax (516) 279-4483 • www.ellman.com

RF ACCESSORIES

Cables

Bipolar Cable

IEC-JX/B

Reusable

For use with Surgitron® Dual RF™, S5, AcuSect™ and RadioLase® Devices (1 pc)



IEC-BC

Sterile/Single-Use

Minor Surgery Bipolar cable

For use with Surgitron® Dual RF™, S5, AcuSect™ and RadioLase® Devices (5 pcs)



JX

Reusable

For use with EMC Devices (1 pc)



JX/E

Reusable

For use with EMC Device (1 pc)



Quick-Connect™ Cable

For use with Suction Tonsil and Suction Adenoid Probes, Disposable Suction Coagulator, Mono Stan RF Forceps, and Rodin™ Micro Surgery handpiece.



IEC-SP020

Reusable

For use with Surgitron® Dual RF™, S5, AcuSect™ and RadioLase® Devices (1 pc)

SP020

Reusable

For use with EMC Device (1 pc)



Snare Cable

For use with Polyp Snare™

IEC-H7B

Reusable

For use with Surgitron® Dual RF™, S5, AcuSect™ and RadioLase® Devices (1 pc)



H7B

Reusable

For use with EMC Device (1 pc)



Lap/Endo Monopolar Cable

IEC-JC

Reusable

For use with Surgitron® Dual RF™, S5, AcuSect™ and RadioLase® Devices (1 pc)

Monopolar cable for Endoscopic instruments



For orders and/or pricing information, please call ellman® at: (800) 835-5355
(516) 594-3333 • fax (516) 279-4483 • www.ellman.com

RF ACCESSORIES

Surg-e®-Doc™



Surg-e®-Doc™
SDS30

The Surg-e®-Doc™ offers the following advantages:

Extra Monopolar port

The Surg-e®-Doc™ offers the ability to utilize an extra monopolar port. This added versatility eliminates the need to change monopolar electrodes within the same handpiece during a procedure.

Simultaneous activation of Surg-e®-Vac™

The Surg-e®-Doc™ extends the life of the vacuum unit by automatically activating the Surg-e®-Vac™ only when the RF energy source is activated. Moreover, the Surg-e®-Vac™ will continue vacuuming for 5 seconds after the RF energy source is deactivated. This allows for automatic removal of residual plume.

Less Clutter

The Surg-e®-Doc™ eliminates the need for separate footswitch activation of Surg-e®-Vac™. This reduces clutter since the only footswitch needed is that of the RF energy source.

For orders and/or pricing information, please call ellman® at: (800) 835-5355
(516) 594-3333 • fax (516) 279-4483 • www.ellman.com

RF ACCESSORIES


Adapters for Surgitron® EMC

Universal Fingerswitch Adapter

EMCKA

Reusable

For use with EMC Device, handpiece not included
(1 pc)



1-Button Fingerswitch Handpiece with Locking Chuck

XK01

Reusable

For use with EMC Universal Fingerswitch Adapter
(2 pcs)



The **ellman®** Universal Fingerswitch Adapter plugs into the side of your Surgitron® FFPF EMC. The adaptor eliminates the need for a foot pedal. The standard foot-controlled handpiece is still operational.



EMC Universal
Fingerswitch Plug

For orders and/or pricing information, please call **ellman®** at: (800) 835-5355
(516) 594-3333 • fax (516) 279-4483 • www.ellman.com

RF ACCESSORIES

For Surgitron® FFPF EMC

EMC Steri Caps™

EMCH35

Reusable

Steri Cap™ (3 for large dial and 3 for small dial) for use with EMC only, 3-Blue and 3-Black (6 pcs)



Bipolar Converter

EMCH95A

For use with EMC Device (1 pc)



Bipolar Converter Steri Caps™

H96

Caps for activation switch for use with Bipolar Converter (3 pcs)



Wall Bracket

HH

Wall Bracket with hardware.

For use with EMC only

Swivel arm and

swivel plate.

(1 pc)



Steri Caps™ enable you to control in total sterility with autoclavable protective dial covers.

The Bipolar Converter provides the surgeon the convenience of switching the **ellman**® Surgitron® EMC from the Handpiece in monopolar mode to the Bipolar Forceps mode with the simple flick of a switch. This eliminates the need to plug and unplug the cords during surgery.

The Wall Bracket space saving accessory attaches easily to any vertical surface.

The **ellman**® Wall Bracket will support up to 40 lbs. and is ideal for the Surgitron® EMC.

For orders and/or pricing information, please call **ellman**® at: (800) 835-5355
(516) 594-3333 • fax (516) 279-4483 • www.ellman.com

RF ACCESSORIES

Handpieces



3-Button Fingerswitch Handpiece with Locking Chuck

IEC-X3FHPB

Reusable (2 pcs)

Accepts 1/16" and 3/32" Shaft Electrodes

For Surgitron® Dual RF™ and S5 Devices



Sterile 3-Button Fingerswitch Handpiece with Locking Chuck

IEC-3FHPB/D

Sterile/Single-Use (10 pcs)

Accepts 1/16" and 3/32" Shaft Electrodes

For Surgitron® Dual RF™ and S5 Devices



3-Button Fingerswitch Blade Handpiece with Locking Chuck

IEC-XBHP3

Reusable (2 pcs)

Accepts Sterile/Single-Use Blades (assorted surgical blades including #10, and #15) and 1/16" Electrode

For Surgitron® Dual RF™ and S5 Devices



1-Button Fingerswitch Handpiece with Locking Chuck

XK01

Reusable (2 pcs)

Accepts 1/16" and 3/32" Shaft Electrodes

For Surgitron® EMC Device



1-Button Fingerswitch Handpiece with Locking Chuck

IEC-X1FHPA

Reusable (2 pcs)

Accepts 1/16" and 3/32" shaft electrodes

One button fingerswitch handpiece for AcuSect™ and RadioLase® Devices

The **ellman®** handpiece is ergonomically designed to fit comfortably in your hand therefore maximizing productivity by reducing operator fatigue and discomfort.

RF ACCESSORIES

Handpieces



3" Length Handpiece with Patented Locking Chuck

IEC-XLP1

Reusable

For Surgitron® Dual RF™, S5, AcuSect™ and RadioLase® Devices

3" Length, Accepts 1/16" and 3/32" shaft, for Laparoscopy and Micro-Larynx probes (2 pcs)



microInsulated Needle Handpiece with Locking Chuck

IEC-XH136 – For Surgitron® Dual RF™, S5, AcuSect™ and RadioLase® Devices (2 pcs)

XH136 – For Surgitron® EMC Device (2 pcs)

Reusable (1 pc). Smaller orifice for insertion of **ellman®**

microInsulated needles only (Cat. No. H137 & H1375)

Included in Cosmetic Treatment package (see page 51)



Blade Handpiece with Locking Chuck

IEC-H15 – Accepts 1/16" Shaft. For Surgitron® Dual RF™, S5, AcuSect™ and RadioLase® Devices

H15 – Accepts 1/16" Shaft. For Surgitron® EMC Device

IEC-QH15 – Accepts 3/32" Shaft. For Surgitron® Dual RF

QH15 – Accepts 3/32" Shaft. For Surgitron® EMC

Reusable (1 pcs). Foot Pedal activation

For use with assorted surgical blades including #10, and #15



Surgitron® EMC
(Monopolar)



Surgitron® Dual RF™, S5 and
RadioLase® (Monopolar)



Foot-Controlled Handpiece with Locking Chuck

IEC-XHP1 – For Surgitron® Dual RF™, S5, AcuSect™ and RadioLase® Devices

XHP1 – For Surgitron® EMC Device

Reusable (2 pcs). Accepts 1/16" and 3/32" shaft



Rodin® Microsurgery Handpiece

XRHP – Length (3.895") 98.93mm, Width (.460") 11.68mm

Reusable (2 pcs). To be used with IEC Quick-Connect™ Cable

Accepts 1/16" Shaft. For Surgitron® Dual RF™, S5, AcuSect™ and RadioLase® Devices

(See page 11 for Quick-Connect™ Cable, Cat. No. IEC-SP020)

(See page 33 for Rodin™ Microsurgery ACe-Tip™ Electrode Set, Cat. No. RACE7)

RODIN® 1

Rodin® RF Microsurgery Set

One RHP Rodin® Microsurgery Handpiece, one IEC-SP020 Quick-Connect™ Cable, one RACE7 Rodin® Microsurgery ACe-Tip™ Set, and instructions for use.



Separate plugs are used in accordance with the corresponding device, the Surgitron® Dual RF™, S5 RadioLase® and the Surgitron® EMC.

For orders and/or pricing information, please call **ellman®** at: (800) 835-5355
(516) 594-3333 • fax (516) 279-4483 • www.ellman.com

RF ACCESSORIES

Bipolar Forceps



J1 Length 9.5cm, (3.7"), Diameter 0.4mm, Straight Tip

ellman® Bipolar Forceps


provide clean wet-field
coagulation.



J1B Length 9.5cm, (3.7"), Diameter 0.4mm, Bent Tip



J2 Length 9.5cm, (3.7"), Diameter 0.3mm, Curved Flat/Angled Tip



J3 Length 7.8cm, (3"), Diameter 0.4mm, Angled Tip



J4 Length 7.8cm, (3"), Diameter 0.3mm, Straight Tip



J5 Length 9.2cm, (3.6"), Diameter 0.3mm, Bent Tip



J8 Length 16.3cm, (6.4"), Diameter 1.0mm, Straight Tip

RF ACCESSORIES

For orders and/or pricing information, please call ellman® at: (800) 835-5355
(516) 594-3333 • fax (516) 279-4483 • www.ellman.com

RF ACCESSORIES

Bipolar Forceps



J9 Length 16.2cm, (6.4"), Diameter 1.0mm, Straight Tip




J10 Length 16.2cm, (6.4"), Diameter 2.0mm, Straight Tip



J11 Length 12.2cm, (4.8"), Diameter 0.4mm, Straight Tip



J12 Length 10.9cm, (4.3"), Diameter 1.0mm, Straight Tip



J14 Length 10.8cm, (4.3"), Diameter 0.7mm, Straight Tip, Insulated, Inner 1mm Exposed



J19 Length 13.6cm, (5.4"), Diameter 2.0mm, Straight Tip



EMC Bipolar Cord
(see page 8)



Surgitron® Dual RF™, S5,
and RadioLase® Bipolar
Cord (see page 8)

© 2002 Ellman International, Inc. All rights reserved.

For orders and/or pricing information, please call ellman® at: (800) 835-5355
(516) 594-3333 • fax (516) 279-4483 • www.ellman.com

RF ACCESSORIES


Bayonet Bipolar Forceps



- J7 Length 17.3cm, (6.8"), Diameter 1.0mm, Bayonet, Flat Tip
J7U Length 17.3cm, (6.8"), Diameter 1.0mm, Bayonet, 45° Up, Flat Tip
J7D Length 17.3cm, (6.8"), Diameter 1.0mm, Bayonet, 45° Down Flat Tip



- J13 Length 15.6cm, (6.1"), Diameter 1.0mm, Bayonet, non-stick epistaxis, Serrated Tip



- J17 Length 17.9cm, (7"), Diameter 1.0mm, Bayonet, single external tip insulated



- J18 Length 16.9cm, (6.0"), Diameter 0.5mm, Bayonet



- J20 Length 16.6cm, (6.5"), Diameter 2.0mm, Bayonet



- J21 Length 19.7cm, (7.8"), Diameter 1.0mm, Bayonet




ellman® Bayonet Bipolar Forceps enhance access to the surgical environment for increased visibility during radiosurgical procedures.

For orders and/or pricing information, please call ellman® at: (800) 835-5355
(516) 594-3333 • fax (516) 279-4483 • www.ellman.com


RF ACCESSORIES

ACe-Tip™ Non-Stick Forceps


Featuring Advanced Composition Alloy "Non-Stick", Low Temperature




ACBF-010 Length 8.9cm, (3.5"), Diameter 0.5mm, Straight Tip



ACBF-011 Length 8.9cm, (3.5"), Diameter 0.5mm, Bent Tip



ACBF-012 Length 11.4cm, (4.5"), Diameter 0.5mm, Straight Tip



ACBF-013 Length 14.0cm, (5.5"), Diameter 0.5mm, Straight Tip



ACBF-014 Length 15.9cm, (6.25"), Diameter 0.5mm, Straight Tip

ACe-Tip™ Technology Advantages:^{*}

- Minimal adherence of tissue
- Low absorption
- Tissue preservation



□ □ □ □ □ □ □ □ □ □ □

For orders and/or pricing information, please call ellman™ at: (800) 835-5355
(516) 594-3333 • fax (516) 279-4483 • www.ellman.com

□ \$ □ □ □ □ □ □ □ □ □ □ □
□ □ □ □ □ □ □ □ □ □ □ □
□ □ □ □ □ □ □ □ □ □ □ □


RF ACCESSORIES

ACe-Tip™ Non-Stick Forceps


Featuring Advanced Composition Alloy "Non-Stick", Low Temperature



ACBF-015 Length 17.8cm, (7.0"), Diameter 0.7mm, Straight Tip




ACBF-016 Length 17.8cm, (7.0"), Diameter 1.2mm, Straight Tip




ACBF-017 Length 17.8cm, (7.0"), Diameter 1.5mm, Straight Tip

ACe-Tip™ Non-Stick Bipolar
Forceps clinical photos



ACBF-021 Length 16.5cm, (6.5"), Diameter 1.0mm, Bayonet, Serrated Tip





ACBF-022 Length 18.4cm, (7.25"), Diameter 1.0mm, Bayonet Tip

For orders and/or pricing information, please call ellman® at: (800) 835-5355

(516) 594-3333 • fax (516) 279-4483 • www.ellman.com

Electrodes



Proximal Shaft Diameter

**Super Fine** - Used when minimal tissue alteration is desired (Diam. 0.005")**Fine** - Used when less tissue alteration is desired (Diam. 0.007")**Regular** - Used to produce more coagulation (Diam. 0.009")

Proximal shaft diameter is available in 3/32" or 1/16" (2.4mm or 1.6mm) to accommodate standard handpiece connection.

ABOUT ellman® Electrodes:

The insulated shaft of most electrodes can be bent or contoured for all surgical procedures.

(Note: Empire microIncision™ Sterile Needles are not bendable)

Yellow electrode tips are Sterile/Single-Use disposables

Blue electrode tips are re-usable

Orange electrode tips are re-usable

Purple electrode tips are re-usable

White electrode tips are a mixture of re-usable and disposable.

See product catalog for details

INFORMATION Regarding Ordering Electrodes:

Code To Electrode Suffixes:

B- Bendable

D- Disposable (Sterile/Single-Use)

Code To Electrode Prefixes:

T- Package of (2 pcs)

X- Package of (5 pcs) (Not all items available in 5-packs)

Most electrodes are available in 1/16" or 3/32" shaft sizes.

Use "Q" as a prefix when ordering 3/32" shaft size, i.e. QTA8B

Longer electrode shaft sizes may be available.

Custom sizes may be requested. (Minimum quantity of 10 each)

Contact ellman's sales representatives for more details.

RF ELECTRODES

Sterile/Single-Use Electrodes

Medical Electrode Set

DSEP40

Sterile/Single-Use

Assortment of 8 styles,

RF electrodes, available in 1/16" Shaft Only (40 pcs per pack)

Includes 5 of each electrode:

A2D, A8D, B1D, B2D, C3D, D8D, E3D, F1D

Sterile/Single-Use Electrodes

Available in 1/16" Shaft Only (25 pcs per pack)

The **ellman®** Sterile/Single-Use RF electrode eliminates the risk of cross contamination as well as ensuring consistent performance and patient confidence.

	A2D	Fine Wire Electrode 3/8" Length
	A3D/OS	Needle Electrode 3/4" Length
	A8D	Vari-Tip™ Wire Electrode Adjustable cutting wire
	B1D	Round Loop Electrode 1/4" Dia.
	B2D	Round Loop Electrode 3/8" Dia.
	C3D	Diamond Loop Electrode 3/16" Dia.
	D8D	Ball Electrode 2.3 mm
	E3D	Scalpel Electrode 1" Length
	E4D	Scalpel Electrode 1/2" Length
	F1D	Broad Needle Electrode
	133D	Pin Electrode 1/8" Length



For orders and/or pricing information, please call **ellman®** at: (800) 835-5355
(516) 594-3333 • fax (516) 279-4483 • www.ellman.com

RF ELECTRODES

Empire® microIncision™ Needles

Empire® microIncision™ Needle Electrodes

Sterile/Single-Use, available in 1/16" or 3/32" (2 per pack)
Use "Q" as a prefix when ordering 3/32" shaft size



4mm
TEE301
TEE302
TEE303
TEE304
TEE305
TEE306
TEE307
TEE308
TEE309
TEE310
TEE311
TEE312
TEE313
TEE314
TEE315

The ellman® Empire®
microIncision™ Needle
Electrodes deliver focused
High Frequency Radiowave
energy.

Empire® microIncision™
Needle Electrodes are used for
precise cutting and coagulation
during radiosurgical procedures.
Proximal shaft diameter is
available in 3/32" or 1/16"
(2.4mm or 1.6mm) to
accommodate standard
handpiece connection.
Empire® microIncision™
Needle Electrodes are supplied
in individual sterile packages
of two, are Single-Use and
should not be re-sterilized.


RF ELECTRODES

Empire® microIncision™ Needles

Empire® microIncision™ Needle Electrodes

Sterile/Single-Use, available in 1/16" or 3/32" (2 per pack)

Use "Q" as a prefix when ordering 3/32" shaft size



Q 1/16" 3/32" 1/8" 5/32" 3/16" 7/32" 1/4" 9/32" 5/16" 11/32" 3/8" 13/32" 7/16" 9/16" 15/32" 1/2" 17/32" 19/32" 5/8" 11/16" 13/16" 15/16" 17/16" 19/16" 1 1/16" 1 1/8" 1 5/16" 1 7/16" 1 9/16" 1 11/16" 1 13/16" 1 15/16" 1 17/16" 1 19/16" 1 1/2" 1 5/8" 1 11/8" 1 13/8" 1 15/8" 1 17/8" 1 19/8" 1 1/4" 1 1/2" 1 5/4" 1 13/4" 1 15/4" 1 17/4" 1 19/4" 1 1/2" 1 5/4" 1 13/4" 1 15/4" 1 17/4" 1 19/4"



TEE316 Length 10cm, 45° angle shaft



TEE318 Length 15cm, 45° angle shaft



TEE319 Length 18cm, 45° angle shaft



TEE320 Length 9cm, curved

TEE321 Length 11cm, curved


TEE322 Length 14cm, curved

TEE323 Length 16cm, curved

TEE324 Bayonet Length 9cm




TEE325 Length (2"), 5cm, curved shaft



TEE327 Length (3/4"), 1.9cm, straight



TEE326 Length (1"), 2.5cm,
40 degree curved shaft



TEE328 Length (1"), 2.5cm,
30 degree curved shaft

Note:
1 inch=2.54cm



Upper Lid Blepharoplasty



For orders and/or pricing information, please call ellman® at: (800) 835-5355
(516) 594-3333 • fax (516) 279-4483 • www.ellman.com

RF REUSABLE ELECTRODES

Reusable Electrodes

Needle Electrodes*

Reusable

1/16" shaft size or 3/32" (2 pcs)

Use "Q" as a prefix when ordering 3/32" shaft size



TA1B Length 3/8", Super Fine Wire

TA2B Length 3/8", Fine Wire

TA9B Length 3/8", Regular Wire

TA3B Length 3/8", Regular Needle


TF1B Length 1/2", Broad Needle

TA8B Vari-Tip™ Microdissection™ Needle
(comes with additional wire and has
an adjustable depth)

T133B Pin Electrode
1/8" Length

Needle Electrodes are commonly used for incisions and excisions in precise surgical procedures.

Choice of Wire Gauge



Round Loop Electrodes*

Reusable

1/16" shaft size or 3/32" (2 pcs)

Use "Q" as a prefix when ordering 3/32" shaft size



TBBB 5/32" Dia. Round Loop, Fine Wire

TB1B 1/4" Dia. Round Loop, Fine Wire

TB4B 1/4" Dia. Round Loop, Regular Wire

TB9B 19/64" Dia. Round Loop, Fine Wire

TB10B 19/64" Dia. Round Loop, Regular Wire

TB2B 3/8" Dia. Round Loop, Fine Wire

TB5B 3/8" Dia. Round Loop, Regular Wire

TB3B 1/2" Dia. Round Loop, Fine Wire

TB6B 1/2" Dia. Round Loop, Regular Wire

Round Loop Electrodes are commonly used for shaving lesions, contouring tissue and obtaining biopsy specimens.

***** ***** ***** ***** *****

For orders and/or pricing information, please call ellman® at: (800) 835-5355
(516) 594-3333 • fax (516) 279-4483 • www.ellman.com

RF REUSABLE ELECTRODES

Reusable Electrodes

Diamond Loop Electrodes*

Reusable

1/16" shaft size or 3/32" (2 pcs)

Use "Q" as a prefix when ordering 3/32" shaft size



- TC1B** 1/8" Diamond Loop, Fine Wire
TC4B 1/8" Diamond Loop, Regular Wire



- TC3B** 3/16" Diamond Loop, Fine Wire
TC6B 3/16" Diamond Loop, Regular Wire



- TC17B** 1/4" Diamond Loop, Fine Wire
TC18B 1/4" Diamond Loop, Regular Wire



- TC19B** 5/16" Diamond Loop, Fine Wire
TC20B 5/16" Diamond Loop, Regular Wire



- TC21B** 7/16" Diamond Loop, Fine Wire
TC22B 7/16" Diamond Loop, Regular Wire

Triangle Loop Electrodes*

Reusable

1/16" shaft size or 3/32" (2 pcs)

Use "Q" as a prefix when ordering 3/32" shaft size



- TC9B** 5/32" Triangle Loop, Fine Wire
TC10B 5/32" Triangle Loop, Regular Wire



- TC7B** 1/4" Triangle Loop, Fine Wire
TC8B 1/4" Triangle Loop, Regular Wire



- TC11B** 19/64" Triangle Loop, Fine Wire
TC12B 19/64" Triangle Loop, Regular Wire



- TC13B** 3/8" Triangle Loop, Fine Wire
TC14B 3/8" Triangle Loop, Regular Wire



- TC15B** 1/2" Triangle Loop, Fine Wire
TC16B 1/2" Triangle Loop, Regular Wire

Diamond Loop Electrodes

are commonly used for

one-stroke elliptical incisions.

Triangle Loop Electrodes

are commonly used for

shaving lesions and

contouring tissue.

For orders and/or pricing information, please call ellman® at: (800) 835-5355
(516) 594-3333 • fax (516) 279-4483 • www.ellman.com

RF REUSABLE ELECTRODES

Reusable Electrodes

Ball Electrodes*

Reusable

1/16" shaft size or 3/32" (2 pcs)

Use "Q" as a prefix when ordering 3/32" shaft size



TD3B 2mm Ball



TD8B 2.36mm Ball



TD9B 3.17mm Ball



TD10B 4.34mm Ball



TD11B 4.75mm Ball

Ball Electrodes are commonly used for coagulation.

Oval Loop Electrodes*

Reusable

1/16" shaft size or 3/32" (2 pcs)


Use "Q" as a prefix when ordering 3/32" shaft size



TP1B 5/32" Oval Loop, Fine Wire

TP2B 5/32" Oval Loop, Regular Wire

Oval Loop Electrodes are commonly used for shaving lesions and contouring tissue.



TP3B 1/4" Oval Loop, Fine Wire

TP4B 1/4" Oval Loop, Regular Wire



TP5B 19/64" Oval Loop, Fine Wire

TP6B 19/64" Oval Loop, Regular Wire



TP7B 3/8" Oval Loop, Fine Wire

TP8B 3/8" Oval Loop, Regular Wire



TP9B 1/2" Oval Loop, Fine Wire

TP10B 1/2" Oval Loop, Regular Wire

For orders and/or pricing information, please call ellman® at: (800) 835-5355
(516) 594-3333 • fax (516) 279-4483 • www.ellman.com

RF REUSABLE ELECTRODES

Reusable Electrodes

Scalpel Electrode*

Reusable

1/16" shaft size or 3/32" (2 pcs)

Use "Q" as a prefix when ordering 3/32" shaft size



TE1B 3/16" Scalpel, 45° Angle

TE2B 3/8" Scalpel


TE3B 1" Scalpel

Scalpel Electrodes are commonly used for incisions and excision of tissue. The blade shape provides broader coagulation than a needle electrode.

Empire® MicroTip Electrodes*

Reusable

2mm exposed tip (2 pcs, one of each)



TEE227 Empire MicroTip Needles
one .010" wire, one .020" wire
Insulated Coated (clear teflon)
Needle Set; 2mm exposed tips

Empire® MicroTip Needles are tapered and have a pinpoint design which produces micro-smooth cutting and coagulation while creating less tissue alteration.

Insulated Needle Electrodes™*

Reusable, 1/2mm exposed tip



D7 Insulated Coated Needle Set (3 pcs)

D6A Fine Insulated Coated Needle .004" Super Fine, White, (1 pc)

D6B Fine Insulated Coated Needle .007" Fine, Purple, (1 pc)

D6C Fine Insulated Coated Needle .009" Regular, Orange, (1 pc)

Insulated Needle Electrodes are used for the treatment of Telangiectasia and for epilation.

Rhinophyma Electrode Set

EE234 Set of 3


Reusable



EE233 19/64" Triangle Loop, (1 pc)

EE232 1/8" Triangle Loop, (1 pc)

EE231 1/2" Triangle Loop, (1 pc)



Rhinophyma

For orders and/or pricing information, please call ellman® at: (800) 835-5355
(516) 594-3333 • fax (516) 279-4483 • www.ellman.com

RF REUSABLE ELECTRODES

Reusable Electrodes

Javate Endoscopic DCR Electrodes*


R6 Javate Endoscopic DCR Set of 7 Assorted Electrodes



□□□□□□□ □□□□□□□ □□□□□□□ □□□□□□□□□*

Mucotome Electrodes*

H140 Mucotome Electrode Set of 1/16" Shaft
Reusable
Set of 3 electrodes, includes H141, H142 & H143
Used for excision of mucosal strip



□□□□□□□ □□□□□□□ □□□□□□□ □□□□□□□□□*

For orders and/or pricing information, please call ellman® at: (800) 835-5355
(516) 594-3333 • fax (516) 279-4483 • www.ellman.com

RF REUSABLE ELECTRODES*

ACe-Tip™ Advanced Composition Alloy

TNAEE286 HLTZ Tonsillotomy Set of 3 ACe-Tip™ electrodes

Includes one each of: TNAEE284, TNAEE285, TNAEE287



TNAEE284 HLTZ Tonsillotomy, ACe-Tip™,
15mm Loop, .010 wire, 1/16" shaft (2 pcs)



TNAEE285 HLTZ Tonsillotomy, ACe-Tip™,
12mm Loop, .010 wire, 1/16" shaft (2 pcs)



TNAEE287 1cm Needle, 3" Curved Sleeve ACe-Tip™, 1/16" shaft (2 pcs)



TNAA9B 3/8" Length
ACe-Tip™ Regular Wire 1/16" shaft (2 pcs)



TNAD3B 2mm Coagulating Ball
ACe-Tip™ 1/16" shaft (2 pcs)



TNATB3B 3mm Tapered Ball
ACe-Tip™ 1/16" shaft (2 pcs)



TNAD13B 5mm Tapered Ball
ACe-Tip™ 1/16" shaft (2 pcs)



TNAB1B 1/4" Round Loop
ACe-Tip™ Fine Wire 1/16" shaft (2 pcs)



TNAB4B 1/4" Round Loop
ACe-Tip™ Regular Wire 1/16" shaft (2 pcs)



TNAB2B 3/8" Round Loop
ACe-Tip™ Fine Wire 1/16" shaft (2 pcs)



TNAB5B 3/8" Round Loop
ACe-Tip™ Regular Wire 1/16" shaft (2 pcs)



TNAA8B Vari-Tip™
ACe-Tip™ 1/16" shaft (2 pcs)
(comes with additional wire and has
an adjustable depth)



TNAEE305 Empire microIncision™
Needle ACe-Tip™ 1/16" shaft (2 pcs)



TNAEE330 Chedly Puctal Occlusion (CPO) Electrodes
Reusable set of two ACe-Tip™ Electrodes

ACe-Tip™ Technology

Advantages:*

- Minimal adherence of tissue
- Low absorption
- Tissue preservation

The ellman® ACe-Tip™ electrodes are designed to reduce resistance/impedance which causes less undesired lateral heat compared to tungsten steel electrodes.

A 45% reduction in the zone of hyalinization has been demonstrated.□

ACe-Tip™ electrodes glide smoothly through tissues and are excellent for sculpting and tissue shaping.

For orders and/or pricing information, please call ellman® at: (800) 835-5355
(516) 594-3333 • fax (516) 279-4483 • www.ellman.com

□ \$ 11,000.00*
\$ 11,000.00+ 100000 100-100
10000" 10 00 0121
□ 1000000 10*

RF REUSABLE ELECTRODES

Electrodes Sets

Note: (See pg. 33 for Disposable Sets)

Rodin® Microsurgery ACe-Tip™ Electrode Set

RACE7

Reusable

1.5" shaft length, 1/16" dia., Bendable Assortment (7 pcs)

Includes one each of:

TREE305, TRA2B, TRA8B, TRD3B, TRBBB, TRC3B, TRB1B



*Standard Electrode Set


S10B

Reusable

1/16" shaft size or 3/32" Use "Q" as a prefix when ordering 3/32" shaft size

Assortment (7 pcs) Includes one each of: C7B, B1B, C3B, A3B, D8B, F1B, A8B

XXXXXXXXXXXX
XXXX6 XXXXX3 XXXXX2&XXXX3 XXXX6 XXXX6 XXXX6 16 16 16 8 6 6 XXXX6 "XXXX6 XXXXX3 XXXXX2&XXXX3 XXXX6 XXXX6 XXXX6 16 16 16 8 6 6 XXXX6 8 6 6 "



Assorted Needle Electrode Set

XA (Formerly XAFB)

Includes one each of: A1, A2, A9, A3, F1 (5 pcs)



Assorted Round Loop Electrode Set

XB (Formerly XBFB)

Fine Wire

Includes one each of: B3, B2, B9, B1, BB

XBRB

Regular Wire

Includes one each of: B6, B5, B10, B4, BB (5 pcs)



For orders and/or pricing information, please call ellman® at: (800) 835-5355

(516) 594-3333 • fax (516) 279-4483 • www.ellman.com

RF REUSABLE ELECTRODES

Electrodes Sets

Assorted Diamond Loop Electrode Set

XC (Formerly XCFB)

Fine Wire

Includes one each of: C21, C19, C17, C3, C1

XCRB Regular Wire

Includes one each of: C22, C20, C18, C6, C4
(5 pcs)



Assorted Ball Electrode Set

XD

Includes one each of: D11, D10, D9, D8, D3
(5 pcs)



Assorted Oval Loop Electrode Set

XP (Formerly XPFB)

Fine Wire

Includes one each of: P9, P7, P5, P3, P1

XPRB

Regular Wire

Includes one each of: P10, P8, P6, P4, P2

(5 pcs)



Assorted Triangle Loops Electrode Set

XT (Formerly XTFB)

Fine Wire

Includes one each of: C15, C13, C11, C9, C7

XTRB

Regular Wire

Includes one each of: C16, C14, C12, C10, C8

(5 pcs)



For orders and/or pricing information, please call ellman® at: (800) 835-5355
(516) 594-3333 • fax (516) 279-4483 • www.ellman.com

RF REUSABLE ELECTRODES

Electrodes Accessories

Electrode Protector Stand™

H66

Electrode Protector Stand™

For regular 2" standard electrodes

Holds 12 electrodes, 6 1/16" x 6 3/32"

(1 pc)


H66A

Electrode Protector Stand™

For large 4" long electrodes

Holds 12 electrodes, 6 1/16" x 6 3/32"

(1 pc)



Electrode Block™

H1

Electrode Block™

Holds 18, 1/16" Electrodes

(1 pc)

H1A

Electrode Block™

Holds 9, 1/16" Electrodes and
9, 3/32" Electrodes

(1 pc)

QH1

Electrode Block™

Holds 18, 3/32" Electrodes
(1 pc)



The **ellman®** Electrode

Protector Stand™ protects
your electrodes from damage
and breakage during the
sterilization process.

The **ellman®** Electrode

Protector Stand™ will stand
on any flat surface for easy
identification and access to
each electrode.

EDUCATIONAL MATERIALS

Educational Accessories

Patient Education Brochures



Pellevé®
Patient Ed
CC09107
(50 pcs)
(8 Page)



Cosmetic
Eye Surgery
CC09076
(50 pcs)
(8 Page)



Facial Veins
CC09116
(50 pcs)
(8 Page)



Remove
Unsightly Skin
CC09025
(50 pcs)
(8 Page)




Treatment
of Snoring
CC10021
(50 pcs)
(8 Page)

Pellevé® Office Pull Up Banners (39.4" Wide X 79" Tall)



Pellevé® Pull Up
CC-PU-13044-1



Pellevé® Pull Up
CC-PU-13044-2



Pellevé® Pull Up
CC-PU-13046



Pellevé® Pull Up
CC-PU-13043

Pellevé® Marketing Kits, Pellevé® Before & After Books and Poster



Pellevé® Marking Kit
CC-PMK-13081 (US version)
and
Pellevé® Marking Kit
CC-PMK-13082 (OUS version)



Pellevé® Before & After Book
CC-BA-13065 (US version)
and
Pellevé® Before & After Book
CC-BA-13066 (OUS version)



Pellevé® Poster
(17" Wide X 26" Tall)
CC-PO-13090

For orders and/or pricing information, please call ellman® at: (800) 835-5355
(516) 594-3333 • fax (516) 279-4483 • www.ellman.com

SMOKE EVACUATION SYSTEM

Device/Accessories



Cart

H157

The elegantly designed H157 Cart includes 2 shelves for accessories, a gooseneck flexible arm, and a convenient front handle for micro positioning. Four rugged caster wheels ensure easy mobility from room to room.

(Shown with Pelleve® S5 Generator, Surg-e®-Vac™, and accessories.-Not included)

NOTE: H157 Cart not compatible with AcuSect, RadioLase or EMC generators.



Surg-e®-Vac™

SVD110 Surg-e®-Vac™ (for 110v countries)

SVD220EU Surg-e®-Vac™ (for 220v countries)



SurgiCart™

H156

SurgiCart™ III with extended base, instrument tray with gooseneck arm, adjustable height. (Shown with Surgitron EMC RF Generator and accessories – Not included)



Vac Assist Stand

W133

The Surg-e®-Vac™ Smoke Evacuation System has been designed with a high suction, high flow rate centrifugal action pump. The ultra-quiet motor is used to draw the surgical smoke from the surgical site through the vacuum tubing and into the Quad Filter™ where the surgical smoke is processed by a series of filters.

NOTE: H156 not compatible with S5 generators.

SMOKE EVACUATION SYSTEM

Accessories

Surg-e®-Vac™ Accessories

The accessories below are used only for Surg-e®-Vac™

Surg-e®-Vac™
Quadtration
Filter
SVF1 (3 pcs)
SVXF1 (12 pcs)



Pneumatic
Footswitch
SVFSC
(10-1/2 Ft. Cable)
(1 pc)



For both Surg-e®-Vac™ and Vapor-Vac II™ Accessories

Plastic Hose **VV148** for
Surg-e®-Vac™ or
VV127 for Vapor-Vac II™
(includes **VV121** Adapter)



Speculum Tubing
Reducing Package
VV130 includes
(2 pcs) **VV128**
(2 pcs) **VV129**



Saf-T-Wand
VV140
(10 pcs)



Speculum Tubing
Reducer Fitting
VV129
(2 pcs)



Speculum Tubing
VV128
1/4" Autoclavable
(2 pcs)



External Tube
VV121
(3 pcs)



Vapor-Vac II™ Accessories

These accessories below are used only for Vapor-Vac II™

Charcoal
Filter
VC123
(2 pcs)



External
Filter
VF124
(10 pcs)



Saf-T-Wand
Filter Set
W146
(25 pcs)



Disposable
External Filter
VV139
(25 pcs)



Vapor-Vac™
Wand Clip
VV131
(1 pc)



For orders and/or pricing information, please call ellman® at: (800) 835-5355
(516) 594-3333 • fax (516) 279-4483 • www.ellman.com

MULTI-SPECIALTY ACCESSORIES

Accessories

Medi-Frig™

Portable Cryosurgery Skin Refrigerant

H11B

1 280g. can, 10 Straw Applicators



XH11B

6 cans, 60 Straw Applicators

H11BCS

25 cans, 250 Straw Applicators

Medi-Frig™ is a portable,

hand-held, cryosurgical

spray unit containing

tetrafluoroethane (-54°C).

The ellman™ cryosurgical spray is safe, storable, easily transportable and cost-effective compared to traditional liquid nitrogen. An added advantage is that it is nonflammable and can be used in conjunction with radiosurgical, electrosurgical and laser procedures.


Wipe-Its™

H19

Pressure sensitive sandpaper pads used to clean radiosurgical electrode tips (4 pcs)

Wipe-Its™ are easy to use,

for effortless removal of tissue from radiosurgery tips.



For orders and/or pricing information, please call ellman® at: (800) 835-5355
(516) 594-3333 • fax (516) 279-4483 • www.ellman.com

MULTI-SPECIALTY ACCESSORIES

Otorhinolaryngology



3" Length Handpiece with Patented Locking Chuck

IEC-XLP1

Reusable

For Surgitron® Dual RF™, S5, AcuSect™ and RadioLase® Devices
3" Length, Accepts 1/16" and 3/32" shaft,
for Laparoscopy and Micro-Larynx probes (2 pcs)

XLP1

Reusable

For the EMC Device
3" Length, Accepts 1/16" and 3/32" shaft,
for Laparoscopy and Micro-Larynx probes (2 pcs)





Micro-Larynx™ RF Probes

IEC-EE275A

Reusable


For Surgitron® Dual RF™, S5, AcuSect™ and RadioLase® Devices
1/16" shaft, 240mm long, (3 pcs); includes IEC-LP1, 3 inch Handpiece



Loop

EE276


1/16" shaft, 240mm long
(1 pc)



Needle

EE277

1/16" shaft, 240mm long
(1 pc)



Ball

EE278A

1/16" shaft, 240mm long
(1 pc)

EE275A for EMC Device

1/16" shaft, 240mm long (3pcs); includes LP1, 3 inch handpiece

The Patented Micro-Larynx™ Probe allows quick and efficient coagulation, and can be used for Granulomas, Papillomatosis, Cyst Biopsies, Vocal Cord Polyps, Posterior Cordectomy, Glottic and Sub-glottic Angiomas, Laryngeal Tumor and Laryngeal Biopsies.

For orders and/or pricing information, please call ellman® at: (800) 835-5355
(516) 594-3333 • fax (516) 279-4483 • www.ellman.com

MULTI-SPECIALTY ACCESSORIES

Otorhinolaryngology

OTO-RAM™ Radiofrequency Assisted Myringotomy Set

IEC-EE600

Reusable

For Surgitron® Dual RF™, S5, AcuSect™ and RadioLase® Devices
Two 4" tubes, 1.5mm and 2.0mm and one Myringotomy Handpiece

IEC-H604

Reusable

For Surgitron® Dual RF™, S5, AcuSect™ and RadioLase® Devices
RF Assisted Myringotomy Handpiece (1 pc)

EE600 and H604 and for EMC Devices

The OTO-RAM™ allows for an effective treatment where office or clinical procedures take place in only a few minutes. Costs for OR Surgery and extensive use of antibiotics and general anesthesia are greatly reduced. The OTO-RAM™ is cost effective compared to laser assisted myringotomy.



2mm hole



2mm hole
Tympanic membrane



OTO-Ram
Probe Tip

OTO-RAM™ Probes

EE601

Reusable

1.5mm Orifice
(1 pc)



EE602

Reusable

2.0mm Orifice
(1 pc)



TEE230

Reusable

Insulated 3.8cm and exposed 5mm
(2 pcs)



For orders and/or pricing information, please call ellman® at: (800) 835-5355
(516) 594-3333 • fax (516) 279-4483 • www.ellman.com

MULTI-SPECIALTY ACCESSORIES

Otorhinolaryngology

Clear-Vu® Bayonet Electrodes

Sterile/Single-Use, (13cm Long Shaft), (1 pc)
1/16" shaft size, 3/32" shaft size available, order with a "Q" prefix



The bayonet design of Clear-Vu® Electrodes provides access to the surgical environment for precise tissue interaction while maintaining a 'clear view' during the procedure.

The assorted tip configurations allows for precise surgical treatment.

Bayonet RF Turbinate Volume Reduction Electrode™

Reusable
Bipolar Turbinate Probe and cord (1 pc)
Tips Insulated ~1.9cm and Exposed ~2.4cm

IEC-16D

For Surgitron® Dual RF™, S5, AcuSect™ and RadioLase® Devices

I6/JXE

For EMC Device



The Bayonet Turbinate Electrode and Bipolar Handpiece design includes a pencil grip for easy handling with a bayonet needle for easy intranasal visualization during surgery.

For orders and/or pricing information, please call ellman® at: (800) 835-5355
(516) 594-3333 • fax (516) 279-4483 • www.ellman.com

MULTI-SPECIALTY ACCESSORIES

Otorhinolaryngology

RAUP Electrode™

H81A1


Sterile/Single-Use
RAUP Radiosurgically Assisted
Uvulopalatoplasty
(2 pcs)



Sub-Mucosal VTR (UPP) Electrode

TEE229A


Sterile/Single-Use
For volumetric reduction
of redundant tissue
1cm insulated, 1 cm exposed
(2 pcs)



Tongue Base Reduction Electrode™

TEE605

Sterile/Single-Use,
For volumetric reduction
of redundant tissue
(2 pcs)



Adenoid RF Electrode™

TEE228

Reusable
1/4" Triangle, 4" Length, Teflon Insulated
(2 pcs)




HTLZ Tonsillotomy Electrode™

TEE284


Reusable
.010 Wire, 15mm Loop
(2 pcs)



Uvulopalatoplasty (RAUP)



Submucosal Stiffening



Tongue Base Procedure

MULTI-SPECIALTY ACCESSORIES

Otorhinolaryngology



Suction Adenoid Probes

IEC-H85

Reusable

Set of 3 Probes, Quick-Connect™ Cable, and Cleaning Brush
(See pg. 11 for Quick-Connect™ Cable)



IEC-H86 Small Probe

9mm with Quick-Connect™ Cable, (1 pc)



IEC-H87 Medium Probe

15mm with Quick-Connect™ Cable, (1 pc)




IEC-H88 Large Probe

19mm with Quick-Connect™ Cable, (1 pc)

Cleaning Brush for Suction Probe

H89

(3 pcs)



Suction Tonsillectomy Probe™

and Quick-Connect™ Cable (See pg. 11 for Quick-Connect™ Cable)

IEC-H83/SP

Reusable

For Surgitron® Dual RF™, S5, AcuSect™ and RadioLase® Devices

H83/SP

For EMC Device



Bipolar Tonsil Reduction Probe™ and Quick-Connect™ Cable

IEC-H84

(See pg. 11 for Quick-Connect™ Cable)

Reusable

For Surgitron® Dual RF™, S5, AcuSect™ and RadioLase® Devices

H84

For EMC Device

For orders and/or pricing information, please call ellman® at: (800) 835-5355
(516) 594-3333 • fax (516) 279-4483 • www.ellman.com

MULTI-SPECIALTY ACCESSORIES


Otorhinolaryngology

Suction Coagulator™

HD41

Sterile/Single-Use

For use with Quick-Connect™ Cable, (12 pcs)
(see page 11 for Quick-Connect™ Cable)




Angled Suction Coagulator™

HD42

Sterile/Single-Use

For use with Quick-Connect™ Cable, (12 pcs)
(see page 11 for Quick-Connect™ Cable)



Suction Coagulation Kit

Reusable

IEC-H40

Shipped with three 3" and three 6" tubes
For Surgitron® Dual RF™, S5, AcuSect™ and
Radiolase® 3 Devices
Includes 1/16", 3/32" and 1/8" Orifice sizes

H40

Shipped with three 3" and three 6" tubes
For EMC Device. Includes 1/16", 3/32"
and 1/8" Orifice sizes



Stylus Round Wire

H44 Suction Tube 6" Long, 1/16" Orifice

H45 Suction Tube 6" Long, 3/32" Orifice

H46 Suction Tube 6" Long, 1/8" Orifice

H47 Suction Tube 3-3/4" Long, 1/16" Orifice

H48 Suction Tube 3-3/4" Long, 3/32" Orifice

H49 Suction Tube 3-3/4" Long, 1/8" Orifice

For orders and/or pricing information, please call ellman® at: (800) 835-5355
(516) 594-3333 • fax (516) 279-4483 • www.ellman.com

MULTI-SPECIALTY ACCESSORIES

Otorhinolaryngology

RF Assisted Tonsillotomy (RAT) Ball Electrodes™

EE264 Set of 3

Reusable

TEE261 2.0mm Ball 4" long (2 pcs)



TEE262 2.36mm Ball 4" long (2 pcs)



TEE263 3.17mm Ball 4" long (2 pcs)




Cryptolysis Electrode

TE1022B

Reusable

(2 pcs)



An RF cryptolysis probe comprises a generally elongated insulated structure leading to a distal end that is angled downward and terminates in a narrow tip that is exposed. The exposed tip is inserted into tonsillar crypts and activated.



The Ball Tonsil Probe is an advancement over the traditional tonsillotomy. The RAT procedure permits controlled reduction of the tonsil without excessive thermal damage or blood loss.



Cryptolysis

Polyp Snare™ System for Nasal Polyps

IEC-H301

Reusable

For Surgitron® Dual RF™, S5, AcuSect™ and RadioLase® Devices with Snare Cable and 3 wires

H301


For EMC Device

with Snare Cable and 3 wires

H7A Snare Wires
6 wires

H7B Snare Cable
For EMC Device (1 pc)

IEC-H7B Snare Cable
For Surgitron® Dual RF™, S5, AcuSect™ and RadioLase® Devices (1 pc)



For orders and/or pricing information, please call ellman® at: (800) 835-5355
(516) 594-3333 • fax (516) 279-4483 • www.ellman.com

MULTI-SPECIALTY ACCESSORIES

General Surgery/Proctology

Laparoscopic Handpiece

IEC-XLP1

Reusable (2 per pack)
3" Length Handpiece.
Accepts 1/16" and 3/32" shafts
For Surgitron® Dual RF™, S5, AcuSect™ and
RadioLase® Devices (1 pc)



IEC-LP1

Reusable (1 ea)

Laparoscopic Electrodes™

Reusable
Lap-Trode Individual Electrodes, 37.5cm Length Electrodes, 1/16" Shaft



L3 2mm Ball Lap-Trode



L4 9.5mm Needle Lap-Trode



L5 2.3mm Hook Lap-Trode



L6 2.3mm Blade Lap-Trode



L8 Hook Needle Lap-Trode



Laparoscopic Electrodes™ Set

L7

Reusable

Lap-Trode Set of 5 Assortment, (L3, L4, L5, L6, L7, L8) 37.5cm Length Electrodes, 1/16" Shaft

MULTI-SPECIALTY ACCESSORIES

Gynecology

LLETZ Electrodes Assortment Set

- W4** One of each in 5 styles, reusable (5 pcs)
W4D One of each in 5 styles, Sterile/Single-Use (5 pcs)

LLETZ Individual Electrodes

Available in 1/16" or 3/32"
Use "Q" as a prefix when ordering 3/32" shaft size
Use "D" as a suffix when ordering Sterile/Single-Use Versions of the LLETZ loop electrodes.

- W1** 15mm x 10mm Loop, reusable (1 pc)
W1D 15mm x 10mm Loop, Sterile/Single-Use (5 pcs)



- W2** 20mm x 10mm Loop, reusable (1 pc)
W2D 20mm x 10mm Loop, Sterile/Single-Use (5 pcs)



- W3** 25mm x 10mm Loop, reusable (1 pc)
W3D 25mm x 10mm Loop, Sterile/Single-Use (5 pcs)



- F2** Endo-Cervical Box 8mm x 11mm Box, reusable (1 pc)
F2D Endo-Cervical Box 8mm x 11mm Box, Sterile/Single-Use (5 pcs)



- W18** 5mm Ball, reusable (1 pc)
W18D 5mm Ball, Sterile/Single-Use (5 pcs)



Conization Electrode™

Reusable



- G1** Conization Electrode
G1/4 4" long, Rotated at 360° to provide a cone shaped specimen

Standard Electrode Set

- S10B/4** 4" shaft length, (7pcs) see page 33 for more information

External Genital Lesion RF Electrodes™

Reusable

- W15** Complete Set of Genital Lesion Electrode set of 4, includes W16, W17, W18 and W19 Electrodes

- W16** 90mm Length, Angled Box



- W17** 4mm Length, Box



- W18** 90mm Length, 5mm Ball



- W19** 90mm Length, 3/8" Needle



The **ellman®** LLETZ loop electrodes are used for the diagnosis and treatment of cervical intraepithelial neoplasia (CIN) and dysplasia.

For orders and/or pricing information, please call **ellman®** at: (800) 835-5355
(516) 594-3333 • fax (516) 279-4483 • www.ellman.com

MULTI-SPECIALTY ACCESSORIES

Gynecology

LLETZ Bayonet Electrodes

- W35** LLETZ Bayonet Electrode Set (5 pcs)
Reusable. Use "D" as a suffix when ordering Sterile/Single-Use
- W35D** Sterile/Single-Use (5 pcs)
- W32** Reusable
Bayonet, 25mm Loop (1 pc)
- W32D** Sterile/Single-Use (5 pcs)
- W31** Reusable
Bayonet, 20mm Loop (1 pc)
- W31D** Sterile/Single-Use (5 pcs)
- W30** Reusable
Bayonet, 15mm Loop (1 pc)
- W30D** Sterile/Single-Use (5 pcs)
- W34** Reusable
Bayonette, 5mm Ball (1 pc)
- W34D** Sterile/Single-Use (5 pcs)
- W33** Reusable
Endo-Cervical Bayonet,
8mm Box (1 pc)
- W33D** Sterile/Single-Use (5 pcs)




Vaginal Speculums

- H55** Vaginal Speculum - Large, 5" Length
- H57** Vaginal Speculum - Standard, 4 1/2" Length
- H58** Vaginal Speculum (2 pcs - Large and Standard)



Vaginal Lateral Wall Retractor

- H69**



Needle Extender/Measured Dose Syringe

- H71**




For orders and/or pricing information, please call ellman® at: (800) 835-5355
(516) 594-3333 • fax (516) 279-4483 • www.ellman.com

Plastic Surgery

Javate/Khan RF Endoscopic Forehead Suction Probe™

Reusable



Javate/Khan Endo Forehead Suction Set

IEC-EE260

Reusable

Includes 2 probes, one Spade, one Scalpel with Suction Handpiece and Cable for Surgitron® 4.0 Dual RF™ and S5 Devices

EE260

For EMC Device

The Javate/Khan Endoscopic Forehead Suction Probe™ utilizes a combination of suction and radiofrequency cutting.

Javate/Khan Coated Flat Scalpel Suction Probe

EE257

Reusable

7/16" Tip Diameter (1 pc)



Javate/Khan Coated Spade Suction Probe

EE258

Reusable

7/16" Tip Diameter (1 pc)



For orders and/or pricing information, please call ellman® at: (800) 835-5355
(516) 594-3333 • fax (516) 279-4483 • www.ellman.com

MULTI-SPECIALTY ACCESSORIES

Plastic Surgery

Patented Dr. Constantin Stan Monopolar Forceps



Micro Stan RF™ Forceps

IEC-MJ21

Reusable

Suggested for Blepharoplasty and Facelift procedures

3.5" Length

Includes Quick-Connect™ Cable and Video CD

For use with Surgitron® Dual RF™ and S5 Devices

MJ21

For EMC Device

The **ellman®** Micro, Mini and Mono Stan Forceps simultaneously dissect and coagulate tissues using the 4.0 MHz radiofrequency.



Mini Mono Stan RF™ Forceps

IEC-MJ11/M

Reusable

Monopolar Forceps with Cutting, Coagulation, Irrigation and Suction Port

Suggested for Facelifts

4.875" Length

Includes Quick-Connect™ Cable and Video CD

For use with Surgitron® Dual RF™ and S5 Devices

MJ11/M

For EMC Device



Facelift



Mono Stan RF™ Forceps

IEC-MJ11

Reusable

Monopolar Forceps with Cutting, Coagulation, Irrigation and Suction Port

Suggested for Breast Augmentation Procedures

6.125" Length

Includes Quick-Connect™ Cable and Video CD

For use with Surgitron® Dual RF™ and S5 Devices

MJ11

For EMC Device

The insides of tips on the Mini Mono and Mono Stan Forceps are serrated.

For orders and/or pricing information, please call **ellman®** at: (800) 835-5355
(516) 594-3333 • fax (516) 279-4483 • www.ellman.com

MULTI-SPECIALTY ACCESSORIES

Dermatology

Cosmetic Treatment Package

For treating telangiectasia and epilation of hair follicles

H135A For EMC Device

IEC-H135A For S5, Surgitron® Dual RF™, AcuSect™ and RadioLase® Devices

Package Includes: One MicroInsulated Needle Handpiece, 50 MicroInsulated Sterile Needles and 50 Patient Education Brochures


XH136 MicroInsulated Needle Handpiece™ For EMC Device (2 pcs)

IEC-XH136 MicroInsulated Needle Handpiece™ For S5, Surgitron® Dual RF™, AcuSect and RadioLase Devices (2 pcs)

H137 MicroInsulated Sterile Needles (50 pcs)

H1375 MicroInsulated Sterile Needles bulk pkg (250 pcs)

CC09116 Patient Education Brochure: FAQ's Facial Veins (50 pcs)



Dermatology Electrode Assortment Set

DDERM40

Assorted, Sterile/Single-Use Electrodes, 1/16" shaft (4 electrode styles)



Insulated Needles

D41 27 Gauge, 0.014" Diameter Needle

D42 30 Gauge, 0.012" Diameter Needle

D43 Insulated Needle Set, 27-Gauge and 30-Gauge (2 pcs, one of each)



Rhinophyma Electrode Set see page 30

Corneal Shields™ see page 55

Chedly Eyelid Retractor see page 55

For orders and/or pricing information, please call ellman® at: (800) 835-5355
(516) 594-3333 • fax (516) 279-4483 • www.ellman.com

The **ellman® Cosmetic Treatment Package** is an effective, comfortable way to remove small unsightly "spider veins" with minimal recovery time, minimal scarring and immediate improvement in appearance.

With the advent of **ellman®** radiowave sclerotherapy, many small facial veins can now be treated without surgery or injections.¹



Telangiectasia Before



Telangiectasia After

1. Ellman, G. Radiowave Sclerotherapy for Telangiectasia. *American Journal of Dermatological Psychiatry*. 2000; 1(1): 1-4.

MULTI-SPECIALTY ACCESSORIES

Pellevé® Wrinkle Reduction System

GlideSafe® Handpieces



P275HPDU GlideSafe® 7.5mm Wrinkle Treatment Handpiece ^{*}
(3 pcs) and Cable (1 pc); Fixed-Life Reusable



P210HPDU GlideSafe® 10mm Wrinkle Treatment Handpiece ^{*}
(3 pcs) and Cable (1 pc); Fixed-Life Reusable



P215HPDU GlideSafe® 15mm Wrinkle Treatment Handpiece ^{*}
(3 pcs) and Cable (1 pc); Fixed-Life Reusable



P220HPDU GlideSafe® 20mm Wrinkle Treatment Handpiece ^{*}
(3 pcs) and Cable (1 pc); Fixed-Life Reusable

□
□ □

□ □

P2CBD Pellevé® Quick-Connect Cable, 2m.
For use with GlideSafe® Handpieces (6 pcs)



DHPG8 Pellevé® Treatment Gel, 8 oz



Temp 62 Fluke 62 Mini Infrared Thermometer

*GlideSafe® handpieces Patent Pending;
only compatible with Pellevé® S5 or
Surgeon® S5 energy sources.

For orders and/or pricing information, please call ellman® at: (800) 835-5355
(516) 594-3333 • fax (516) 279-4483 • www.ellman.com


MULTI-SPECIALTY ACCESSORIES

Ophthalmology/Oculoplastics

Goisis RF Elevator ACe-Tip™

GE1

Reusable
(1 pc)



Wedge Resection Electrode Set™

EE255

Reusable
(set of three)

Wedge Resection Electrode

TEE237 1/4" tip (2 pcs)



TEE236 3/16" tip (2 pcs)



TEE235 1/8" tip (2 pcs)



Beigi Minimally Invasive Entropion Electrode

TEE603

Sterile/Single-Use
(2 pcs)



3cm

Ophthalmic Micro-Dissection Electrodes™



TA10/C1

Reusable. 2" Electrodes
1.5mm of active .004 Wire (2 pcs)

TA10/C2

Reusable. 2" Electrodes
3.0mm of active .004 Wire (2 pcs)

TA10/C3

Reusable. 2" Electrodes
5.0mm of active .004 Wire (2 pcs)

A10

Reusable

Set of two (1.5mm and 3.0mm, .004 Wire - one of each)

Note: Add/OP after the Catalog Number to order 1/2" shorter,
Ophthalmic Micro-Dissection Electrodes™ (Not available for all codes)

For orders and/or pricing information, please call ellman® at: (800) 835-5355
(516) 594-3333 • fax (516) 279-4483 • www.ellman.com

Goisis RF Elevator ACe-Tip™

Product Benefits:

- Aesthetic midface lift easily done under local anaesthesia.
- Easier and faster orbital decompression in Basedow disease.
- Easier and faster treatment of orbital and facial traumas.
- Versatility – evisceration of the eye with a single instrument.
- Conjunctival peritomy
- Cornea excision
- Globe content removal and full-thickness sclerotomies.

MULTI-SPECIALTY ACCESSORIES

Ophthalmology/Oculoplastics

Corneal Shields™

H781

Autoclavable
Assorted Package
one pair each size,
small, medium, large and
1 suction handle



H75A Small - Red (1 handle, 2 pair)

H76A Medium - Green (1 handle, 2 pair)

H77A Large - Blue (1 handle, 2 pair)

H78A Suction Handle (2 pcs)

H78B Suction Handle (25 pcs)



Chedly Eyelid Retractor

H21A

Reusable
(1 pc)



The **ellman® Corneal Shield™** tightens the skin to allow for more precise incisions. It is an excellent aid when performing periorbital surgery and aids in the protection from cautery, desiccation, irritation or trauma, during and after surgery of eyelid and orbital regions.

Chedly Punctal Occlusion (CPO) Electrodes

TNAEE330

Reusable set of two ACe-Tip™ Electrodes



The Chedly Eyelid Retractor allows access for oculoplastic procedures i.e. upper and lower blepharoplasty, ptosis, trichiasis and minor eyelid procedures.

Baker Brow Probe

EE279

Reusable
(2 pcs)



For orders and/or pricing information, please call **ellman®** at: (800) 835-5355
(516) 594-3333 • fax (516) 279-4483 • www.ellman.com

MULTI-SPECIALTY ACCESSORIES

Family Practice/Podiatry



Digit Tourniquets

H97 ellman® Cots – Small (12 pcs)

H98 ellman® Cots – Medium (12 pcs)

H99 ellman® Cots – Large (12 pcs)

ECOTS - 4 each small, medium and large ellman® Cots.



Ingrown Toenail

The ellman® radiosurgical nail matrixectomy technique for Ingrown Toenail allows the matrix to be destroyed electrically with a teflon coated electrode, which protects the overlying eponychium.

Insulated Matrixectomy Electrodes™*

H9 Reusable

Set of 3 electrodes

1/16" or 3/32" shaft

Matrixectomy Electrode Set

Includes one each of: H10D, H10E, H10A

H10D 2mm Matrix Electrodes
Reusable (2 pc)

H10E 4mm Matrix Spicule Electrodes
Reusable (2 pc)

H10A Matrix Spicule (Ball) Electrode
Reusable (1 pc)

1 2 3 4 5 6 7 8 9 10 11 12 13 14 15 16 17 18 19 20 *

For orders and/or pricing information, please call ellman® at: (800) 835-5355
(516) 594-3333 • fax (516) 279-4483 • www.ellman.com

INDEX

INDEX	5
(")	
(")	
Laser Office Pull Up Banners	35
Patient Educational Brochures	35
Pellevé® Before and After Book	35
Pellevé® Marketing Kits	35
Pellevé® Posters	35
D ")	36
Cart	36
SugiCart™	36
Surg-e®-Vac™	36
Vac Assist Stand	36
D ")	37
Surg-e®-Vac™ Accessories for both Surg-e®-Vac™ and Vapor-Vac II™ Accessories	37
D ")	38
Medi-Frig™	38
Wipe-Its™	38
D ") - D ")	39
D ") , D ")	40
Bipolar Converter Steri Caps™	14
EMC Steri Caps™	14
Patented 1-Button Fingerswitch Handpiece with Locking Chuck	13
Patented Bipolar Converter	14
Universal Fingerswitch Adapter	13
Wall Bracket	14
3 #4	41
3 #4	42
D ")	43
Bipolar Cable	11
Lap/Endo Monopolar Cable	11
Snare Cable	11
Quick-Connect™ Cable	11
D ")	44
Cosmetic Treatment Package	51
Dermatology Electrode Assortment Set	51
Insulated Needles	51
D ")	45
(")	46
Empire® microIncision™ Needle Electrodes	24-25
Medical Electrode Set	23
9 ")	26
ACe-Tip™ Advanced Composition Alloy	31
Ball Electrodes	28
Chedly Puctal Occlusion (CPO) Electrodes	31
Diamond Loop Electrodes	27
Empire® MicroTip Electrodes	29
Insulated Needle Electrodes™	29
Javate Endoscopic DCR Electrodes	30
Needle Electrodes	26
Oval Loop Electrodes	28
Radiosurgical Mucotome Electrodes	30
Rhinophyma Electrode Set	29
Round Loop Electrodes	26
&	
Scalpel Electrode	29
Triangle Loop Electrodes	27
9 ")	32
Assorted Ball Electrode Set	33
Assorted Diamond Loop Electrode Set	33
Assorted Oval Loop Electrode Set	33
Assorted Round Loop Electrode Set	32
Assorted Triangle Loops Electrode Set	33
Rodin® Microsurgery ACe-Tip™ Electrode Set	32
Standard Electrode Set	32
(")	33
AcuSect®	8
Pellevé® S5	5
Surgeon® Dual RF™ S5	6
Surgeon® Dual RF™	7
Surgeon® Dual/EMC	7
Surgeon® FFPF EMC	8
0 #	34
ellman® Digit Tourniquets	55
Insulated Matrixectomy Electrodes™	55
0	35
Dual Frequency Footswitch and Cable Pneumatic Footswitch	10
RadioLase I Footswitch and Cable	10
RadioLase 3 Monopolar/Bipolar Footswitch	10
: ")	46
Laparoscopic Electrodes™	46
Laparoscopic Handpiece	46
: 4	47
Conization Electrode™	47
External Genital Lesion RF Electrodes™	47
LLETZ Bayonet Electrodes	48
LLETZ Electrodes Assortment Set	47
LLETZ Individual Electrodes	47
Monsel's Solution	48
Needle Extender/Measured Dose Syringe	48
Vaginal Lateral Wall Retractor	48
Vaginal Speculums	48
0 #	49
3" Length Handpiece with Patented Locking Chuck	16
Blade Handpiece with Locking Chuck	16
Foot-Controlled Handpiece with Locking Chuck	16
microInsulated Needle Handpiece with Locking Chuck	16
1-Button Fingerswitch Handpiece with Locking Chuck	15
3-Button Fingerswitch Handpiece with Locking Chuck	15
Blade Handpiece with Locking Chuck	15
3-Button Fingerswitch Handpiece with Locking Chuck	15
Sterile 3-Button Fingerswitch Handpiece with Locking Chuck	15
Rodin® Microsurgery Handpiece	16
0 #	50
Javate/Khan Coated Flat Scalpel	49
Javate/Khan Coated Spade	49
Javate/Khan Endo Forehead Suction Set	49
Micro Stan RF™ Forceps	50
Mini Mono Stan RF™ Forceps	50
Mono Stan RF™ Forceps	50
0	51
EMC Neutral Plate	9
IEC Neutral Plate	9
IEC Disposable Neutral Plate-Small Size	9
IEC Disposable Neutral Plate-Standard Size	9
0 ") < = ") / ")	52
GlideSafe® Handpieces	52
Pellevé® Quick-Connect Cable	52
Pellevé® Products & Accessories	52
> , ")	53
Baker Brow Probe	55
Beigi Minimally Invasive Entropion Electrode	54
Chedly Eyelid Retractor	55
Chedly Puctal Occlusion (CPO) Electrodes	55
Corneal Shields™	55
Gosis RF Elevator ACe-Tip™	54
Ophthalmic Micro-Dissection Electrodes™	54
Wedge Resection Electrode	54
Wedge Resection Electrode Set™	54
> ")	55
3" Length Handpiece with Locking Chuck	39
Adenoid RF Electrode™	42
Bayonet RF Turbinate Volume Reduction Electrode™	41
Cleaning Brush for Suction Probe	43
Clear-Vu® Bayonet Electrodes	41
Cryptolysis Electrode	45
HLTZ Tonsillectomy Electrode™	42
Micro-Larynx™ RF Probes	39
OTO-RAM™ Radiofrequency Assisted Myringotomy Set	40
OTO-RAM™ Probes	40
RAUP Electrode™	42
Angled Suction Coagulator™	44
Bipolar Tonsil Reduction Probe™ and Quick-Connect™ Cable	43
Myringotomy Electrode™	40
Suction Adenoid Probes	43
Suction Coagulator™	44
Suction Coagulation Kit	44
Suction Tonsillectomy Probe™ and Quick-Connect™ Cable	43
Tongue Base Reduction Electrode™	42
Polyp Snare™ System for Nasal Polyps	45
Quick-Connect™ Cable	39
RF Assisted Tonsillectomy (RAT) Ball Electrodes™	45
Sub-Mucosal VTR (UPP) Electrode	42
0 #	56
Javate/Khan Coated Flat Scalpel	49
Javate/Khan Coated Spade	49
Javate/Khan Endo Forehead Suction Set	49
Micro Stan RF™ Forceps	50
Mini Mono Stan RF™ Forceps	50
Mono Stan RF™ Forceps	50
0	57

For orders and/or pricing information, please call ellman® at: (800) 835-5355
 (516) 594-3333 • fax (516) 279-4483 • www.ellman.com

NOTES

For orders and/or pricing information, please call ellman® at: (800) 835-5355
(516) 594-3333 • fax (516) 279-4483 • www.ellman.com

NOTES

For orders and/or pricing information, please call ellman® at: (800) 835-5355
(516) 594-3333 • fax (516) 279-4483 • www.ellman.com



Ellman International, Inc.
1000 Corporate Park Drive • Suite 100 •
Brentwood, CA 94513 • 925.467.1100 •
www.ellman.com

CC13156 Rev B

Copyright 2013 ellman® International, Inc. All Rights Reserved.

Note: Do not copy or distribute without written authorization from Ellman International Inc.

