



**veterled**   
*animals' friend.*

*vecar / veterinary tecar system.*



## vecar, tecar therapy.

**Vecar** was born from the sports-medical experience and high-level veterinary medicine, for non-invasive orthopedic rehabilitation treatments, through the generation of endogenous energy concentrating it where you really want to apply it. With the use of different application techniques (resistive or capacitive) and specific configurations you can get the kind of action that you want on the tissue at both superficial and deep level as the endogenous energy produced does not undergo variations due to the natural barrier of the skin.

Among the effects we have:

- Intra-and extracellular Biostimulation
- Hyperthermia
- Hyperemia
- Drainage post-traumatic
- Powerful relaxant action
- Effective analgesic action

The most advanced technology combined with the knowledge of medical-veterinary, allows the direct treatment of target tissues, previously difficult to access, obtaining the best results without the use of drugs and in a completely non-invasive way.

01 configuration 0°



02 configuration 45°



03 configuration 90°



### characteristics

Max Power  
**150 W**  
 Working Frequencies  
**500 - 1000 kHz**  
 Operative mode  
**Resistive/Capacitive**  
 Power supply  
**90-240 Vac / 50-60 Hz**  
 Electrical input power  
**300 VA**




### novità

- 7 Inches Touchscreen Display
- User Friendly Interface Software
- Updates via USB ports
- Control of the device via the display touchscreen, encoder and using the keyboard
- Selectable Protocols
- Customizable programs
- Security System Security slot
- 6 selectable languages
- Handpiece with three work configurations and three interchangeable heads with different diameters



Size (mm)  
 width 315 / height 170 / depth 390  
 weight 5 kg

**vecar,**  
*animals' friend.*





ASKMED